

SANITARY ANGLE BLEED AND SAMPLE VALVES (Safety Cap)

BENEFITS

- Minimal Dead Leg
- Fully Autoclavable / Sanitary design
- Indicating Stem for Open / Close
- Non Rotating Stem
- Quick and easy assembly / disassembly per your maintenance schedule
- Fully customizable to customer specifications

TECHNICAL INFORMATION

- Controlled sulphur stainless steel
- Available in Alloy 22
- Variety of surface finishes available
- Selection of "O" rings available:
EPDM - Viton - Silicone (Clear/Red)

- PTFE-TFM Seat
- Variety of chemical treatments available
- Variety of Valve sizes:
 - 1/8" valve = 0.104" Inlet Bore size
 - 1/4" valve = 0.187" Inlet Bore size
 - 1/2" valve = 0.370" Inlet Bore size
- MTR'S standard for full material traceability (Heat numbers/codes)



www.fluidlinetech.com

