


Approximate dimensions of STANDARD, FLANGED and ENVELOPE style gaskets.


Approximate dimensions of STANDARD Tri-Clamp style gaskets

DIMENSIONS	ID	OD	
1/2"	3/8"	1"	For 1/2" and 3/4" Approx. thickness: .079"
3/4"	5/8"	1"	
1"	7/8"	2"	
1 1/2"	1 3/8"	2"	
2"	1 7/8"	2 1/2"	
2 1/2"	2 3/8"	3"	
3"	2 7/8"	3 1/2"	
4"	3 7/8"	4 1/2"	
5"	4 3/4"	5 1/2"	
6"	5 3/4"	6 1/2"	



Approximate dimensions of STANDARD FLANGED Tri-Clamp style gaskets

DIMENSIONS	ID	OD	
1"	7/8"	2"	
1 1/2"	1 3/8"	2"	
2"	1 7/8"	2 1/2"	
2 1/2"	2 3/8"	3"	
3"	2 7/8"	3 1/2"	
4"	3 7/8"	4 1/2"	
5"	4 3/4"	5 1/2"	
6"	5 3/4"	6 1/2"	
8"	7 3/4"	8 1/2"	
10"	9 11/16"	10 1/2"	
12"	11 11/16"	12 1/2"	



Approximate dimensions of TEFLON ENVELOPE Tri-Clamp style gaskets

DIMENSIONS	ID	OD	
1/2"	3/8"	7/8"	
3/4"	5/8"	7/8"	
1"	7/8"	2"	
1 1/2"	1 3/8"	2"	
2"	1 7/8"	2 1/2"	
2 1/2"	2 3/8"	3"	
3"	2 7/8"	3 1/2"	
4"	3 7/8"	4 1/2"	
5"	4 3/4"	5 1/2"	
6"	5 3/4"	6 1/2"	
8"	7 3/4"	8 1/2"	
10"	9 3/4"	10 1/2"	
12"	11 3/4"	12 1/2"	

