


# Sanitary hangers for tube, hose & pipe



## Patented Rounded Hygienic Hangers

- Patented Design to latest BPE standard
- 316 polished stainless housings, silicone or PTFE inserts
- Easily cleaned with minimized areas for contaminants to be trapped
- Compatible with steam line temperatures to 400° F (204°C)
- Flexible inserts allows required slope for draining
- Multiple designs to mounting from ceiling, floor, wall or be stacked



## Smooth Bore Block Style Supports

- Santoprene, Polypropylene and PTFE inserts for broad thermal range
- Available as anchors or thermal expansion guides
- Multiple mounting configurations for welding, anchoring hanging, stacking... etc



## Slim Line Hygienic Design

- Captive fasteners – NO LOOSE BOLTS
- Fast and easy assembly
- 304 polished stainless housings
- Easily cleaned, no hinges or crevasses



## Hex Hangers – Sanitary Supports

- Tube sizes to 6"
- EPDM grommets and Polypropylene sleeves
- Consult factory for weld couplings and mounting rods



## Hygienic Clamp Unions

- Single pin heavy duty up to 6"
- High pressure 2 bolt design
- Three piece design

7 Wm. Demarest Place  
Waldwick, NJ 07463

 **STAUFF**<sup>®</sup>  
PROCESS PRODUCTS

201.444.7800  
sales@stauffusa.com