



Top-Flo[®] C-Series Pump Dolly



Stainless Steel Flow Control Equipment for the
Food, Beverage, Dairy, Cosmetics, Biotechnology,
Pharmaceutical and Electronics Processing Industries

Top-Flo® C-Series Pump Dolly

TOP-FLO® C-Series pump dolly has an emphasis on your convenience. Its mobility will allow you to find many functions in your manufacturing process.

Standard features are as follows:

- Universal dolly for C100, C114, C216, C218 & C328 pumps with motors through 10 horsepower
- Full 304 stainless steel construction with 2 fixed wheels
- Bead blast finish
- 8" semi-pneumatic wheels
- Adjustable pivot point axle
- Handle mounted cord wrap

Optional features are as follows:

- Finishes:
 - Mechanical polish
 - Electropolish
- Full range of control options:
 - VFD
 - Start/stop controls
 - Network controls
- Wheels:
 - Non-marking white wheels
 - Full pneumatic
 - Solid
- Special customer design requirements

Advantages:

- Designed for your convenience
- Can be designed & manufactured to your specifications
- Convenient roll-away design
- Mobility allows you to find many functions in your manufacturing process
- Durable (full 304 stainless steel construction)



. . . the source for
all your product needs
in the process industries

Contact your Top Line Representative for Assistance
Office: P.O. Box 264 · Bradford, PA 16701

Plant: 21 Valley Hunt Drive · Lewis Run, PA 16738

Phone: 800-458-6095

814-362-4626

Fax: 814-362-4453

www.toplineonline.com